
COLLOIDAL IRON - MUCOPOLYSACCHARIDES

PURPOSE: To demonstrate carboxylated and sulfated mucopolysaccharides and glycoproteins.

PRINCIPLE: Colloidal ferric ions are, at low pH, absorbed principally by carboxylated and sulfated mucosubstances. The excess reagent is washed out and the classic Prussian blue reaction (iron stain) is used to demonstrate iron bound to the tissue.

CONTROL: Use mucin control, small intestine

FIXATION: 10% formalin, avoid chromate fixatives.

TECHNIQUE: Cut paraffin sections at 4m.

EQUIPMENT: Hot plate, acid-clean glassware: flasks, cylinders, and coplin jars, filter paper.

REAGENTS:

Müller Colloidal Iron

Stock Solution:

Boiling distilled water 250.0 ml
29% ferric chloride-fresh 4.4 ml

Keep the water boiling and stirring when adding ferric chloride. Stir until solution becomes dark red. Remove from heat, cool. Store in clean brown bottle, label with date and initial. Solution is stable for 6 months.

Working Solution:

Stock colloidal 10.0 ml
Distilled water 18.0 ml
Acetic acid 12.0 ml

Prepare fresh, discard after use.

CAUTION: Avoid contact and inhalation.

12% Acetic Acid:

Glacial acetic acid 24.0 ml
Distilled water 176.0 ml

CAUTION: Corrosive acid.

Ferrocyanide-Hydrochloric Acid Working Solution (Iron Stain)

5% Potassium ferrocyanide 25.0 ml
5% Hydrochloric acid 25.0 ml

Mix right before use, use acid clean glassware. Discard after use.

CAUTION: Corrosive acid.

Nuclear-Fast Red:

See Iron Stain

CARBOHYDRATES

COLLOIDAL IRON - MUCOPOLYSACCHARIDES

Page: 2 of 2

SAFETY: Wear gloves, goggles and lab coat. Boil solution in fume hood. Add acid to water. Avoid contact and inhalation.

Acetic acid: Irritating to eyes, skin and respiratory system, target organ effects on respiratory system. Corrosive.

Ferric chloride: corrosive to metal. Causes GI problems on ingestion.

PROCEDURE:

1. Deparaffinize and hydrate sections to distilled water.
2. Rinse slides in 12% acetic acid.
3. Working colloidal iron solution for 1 hour.
4. Rinse in 12% acetic acid, 3 changes, 3 minutes each.
5. Ferrocyanide-hydrochloric acid solution for 20 minutes.
6. Wash in running tap water for 5 minutes.
7. Counterstain with Nuclear-fast red for 5 minutes.
8. Wash in running tap water.
9. Dehydrate, clear and mount.

RESULTS:

Acid mucopolysaccharides and sialomucins	deep blue
Nuclei	pink-red
Cytoplasm	pink

REFERENCE:

F.CARSON, HISTOTECHNOLOGY: A SELF-INSTRUCTIONAL TEXT, PP 127-129, 1991, ASCP PRESS

Crookham,J, Dapson,R, Hazardous Chemicals in the Histopathology Laboratory, 2nd ED, 1991, Anatech

Prepared: _____ By: _____

Approved: _____ By: _____

Downloaded from WebPath: Internet Pathology Laboratory
<http://www-medlib.med.utah.edu/WebPath/webpath.html>

PROCEDURE CARD
COLLOIDAL IRON - MUCOPOLYSACCHARIDES

CONTROL: Use mucin control, small intestine

PROCEDURE:

1. Deparaffinize and hydrate sections to distilled water.
2. Rinse slides in 12% acetic acid.
3. Working colloidal iron solution for 1 hour.
4. Rinse in 12% acetic acid, 3 changes, 3 minutes each.
5. Ferrocyanide-hydrochloric acid solution for 20 minutes.
6. Wash in running tap water for 5 minutes.
7. Counterstain with Nuclear-fast red for 5 minutes.
8. Wash in running tap water.
9. Dehydrate, clear and mount.

RESULTS:

Acid mucopolysaccharides and sialomucins	deep blue
Nuclei	pink-red
Cytoplasm	pink

REAGENTS:

Muller Colloidal Iron

Stock Solution:

Boiling distilled water	250.0 ml
29% ferric chloride-fresh	4.4 ml

Keep the water boiling and stirring when adding ferric chloride. Stir until solution becomes dark red. Remove from heat, cool. Store in clean brown bottle, label with date and initial. Solution is stable for 6 months.

Working Solution:

Stock colloidal	10.0 ml
Distilled water	18.0 ml
Acetic acid	12.0 ml

Prepare fresh, discard after use.

Ferrocyanide-Hydrochloric Acid

Working Solution

(Iron Stain):

5% Potassium ferrocyanide	25.0 ml
5% Hydrochloric acid	25.0 ml

Mix right before use, use acid clean glassware. Discard after use.

CAUTION: Corrosive acid.

12% Acetic Acid:

Glacial acetic acid	24.0 ml
Distilled water	176.0 ml

CAUTION: Corrosive acid.

Nuclear-Fast Red:

See Iron Stain

MULLER COLLOIDAL IRON

Stock Solution:

Boiling distilled water 250.0 ml
29% ferric chloride-fresh 4.4 ml

Keep the water boiling and stirring when adding ferric chloride. Stir until solution becomes dark red. Remove from heat, cool. Store in clean brown bottle, label with date and initial. Solution is stable for 6 months.

Working Solution:

Stock colloidal 10.0 ml
Distilled water 18.0 ml
Acetic acid 12.0 ml

Prepare fresh, discard after use.

CAUTION: Avoid contact and inhalation.

TECH: _____

DATE: _____

EXPIRATION: _____

12% ACETIC ACID:

Glacial acetic acid 24.0 ml
Distilled water 176.0 ml

Mix well, label with initial and date. Solution is stable for 6 months.

CAUTION: Corrosive acid.

TECH: _____

DATE: _____

EXPIRATION: _____