

Subject:

ZENKER'S SOLUTION

PURPOSE: Zenker's is a nuclear fixative and is used as a mordant in some staining procedures.

REAGENTS:

Zenker's Stock Solution:

Mercuric chloride 50.0 gm
Potassium dichromate 25.0 gm
Sodium sulfate 10.0 gm
Distilled water 1000.0 ml

Working Solution:

Zenker's stock 95.0 ml
Acetic acid 5.0 ml

Mix well, stable for 2 years.

Mix just before use, discard following hazardous waste regulations.

CAUTION: CARCINOGENIC, CORROSIVE

SAFETY/PPE: Wear gloves, goggles and lab coat. Work in well ventilated area. Keep hot solution under fume hood.

Mercuric chloride: Severe skin and eye irritant; target organ effects on reproductive, urogenital, respiratory, gastrointestinal and fetal systems via ingestion and inhalation routes.

Potassium dichromate: toxic by inhalation of dust and by ingestion; target organ effects on reproductive and fetal systems. Solid form is corrosive to eyes, skin and mucous membranes.

Acetic acid: Skin and eye irritant; target organ effects on respiratory system by inhalation.

DISPOSAL: Pour used Zenker's into glass bottle, to be disposed of by Safety Services, attach Hazardous waste label. NOTE: "Contains Mercury".

PROCEDURE:

1. The tissues must be treated for mercury pigment.
2. The tissues must be washed well in running water to remove the excess potassium dichromate or the dehydrating alcohols will cause a chrome pigment.

REFERENCE: Carson,F, Histotechnology A Self-Instructional Text, pp19-20, 1st Ed, 1990, ASCP Press

Crookham,J, Dapson,R, Hazardous Chemicals in the Histopathology Laboratory, 2nd ED, 1991, Anatech

Prepared: Sept. 1993 By: Susan Wall

Approved: _____ By: _____

Zenker's Stock Solution:

Mercuric chloride	50.0 gm
Potassium dichromate	25.0 gm
Sodium sulfate	10.0 gm
Distilled water	1000.0 ml

Mix well, stable for 2 years.

Working Solution:

Zenker's stock	95.0 ml
Acetic acid	5.0 ml

Mix just before use. Do not dispose of down the drain, pour into glass bottle labeled "mercury waste".

CAUTION: CARCINOGENIC, CORROSIVE

DATE: _____

TECH: _____

EXPIRATION: _____